

What Does It Measure?

The MoistureMap MM 200 is a device, featuring a capacitance based sensor comparable to the renowned L'Oréal SkinChip®*. The sensor gives graphic information on the **near surface hydration distribution** and the **micro-topography of skin and other tissues** (textiles, plants, etc.).

The Measuring Principle

The sensor measures the **penetration of the electromagnetic field**. On the 18.0 x 12.8 mm silicon chip of the sensor, over 90,000 capacitors are located. Conductive material e.g. water will reflect the signal making the resulting pixel darker while non-conductive material will make the signal penetrate farther inside, and the resulting pixel will be lighter on a scale of 255 grey levels. Rather than absolute moisture figures, the MoistureMap indicates the distribution of hydration on the skin surface. With a special image analysis software, the image can be evaluated in different ways.

Fields of Application

Wherever moisture distribution plays a role, the MoistureMap MM 200 is an ideal imaging addition to the purely quantitative measurements.



- **Efficacy testing** of cosmetics & pharmaceuticals & surfactants
- Sun damage and **illustration** of skin lesions and scars.
- To map the **hair moisture level**

Advantages

- **Flat sensor surface** (no metal rim) for comfortable placement on the skin.
- Quick and easy to handle: **image capture/live image** by a button on the probe handle.
- **Spring** loaded sensor
- **In-vivo** skin measurement and also **in-vitro** applications can be performed.
- **Calibrated probe**, no further calibration by the user.

- **Automatic saving** of the images (.jpg) under study name
- **Video** possible (.avi).
- **Evenness** of the hydration is displayed in 5 different colours and a histogram.
- Additionally **topographic measurements** (profile, corner density, anisotropy index) give interesting aging parameters.
- Selection and comparison of 1-2 **regions of interest, (ROI)**.
- Impressive **“3D”-style images**.
- Easy **export of all images** and results to Excel®.
- Up to six images together with their complete results can be **compared** in one overview.
- The only instrument working side by side with the established **Corneometer®** and **Tewameter®**.

*The MoistureMap is licensed worldwide under the L'Oréal patent for the Skin Chip® (EP 1 438 922 B1). A variety of articles on the measurement principle of the Skin Chip® (same as MoistureMap) has been published.



Technical Data
 Device: Dimensions: 13 x 14.6 x 5 cm, Weight: approx. 1.5 kg, Power supply: external 100-240 VAC, 47-63 Hz, DC 12V/4A, Port: USB 2.0, type B connector
 Probe: Dimensions: length: 16.6 cm, measurement head: 4.3 x 3 cm, Weight: approx. 90 g, Active measurement area: 18,0 x 12.8 mm, Sensor size: 256 x 360 pixel, Sensor resolution: 508 DPI 8Bit/pixel, Measurement principle: relative permittivity; MoistureMap in-vitro Adapter (optional): Dimensions: 23 cm (H) x 8 cm x 8 cm, Weight: 220 g
 Technical changes may be made without prior notice.

MONADERM
 5, rue des Violettes
 98000 MONACO
 ☎ : +377 93 25 26 08
 contact@monaderm.com

